NANOCOLOR® VARIO

Bedienungsanleitung  
User manual  
Mode d'emploi  
Manual del usuario

**NANOCOLOR® VARIO 4**  
**NANOCOLOR® VARIO C2**

**MACHERY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)

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## 1. Technical Data

<b>NANOCOLOR® VARIO 4:</b>	Programmable heating block for chemical-analytical digestions with 2 independent heating units with 12 holes each for test tubes with 16 mm outer diameter
<b>NANOCOLOR® VARIO C2:</b>	Programmable heating block for chemical-analytical digestions with 12 holes for test tubes with 16 mm outer diameter
Display:	Coloured, backlit LCD touchscreen
Operation:	Display menu guidance via touchscreen
Temperature range:	40 – 160 °C (1 °C increments)
Temperatures:	5 pre-programmed temperatures: 70 / 100 / 120 / 148 / 160 °C 7 free memory locations for individual temperature settings
Temperature stability:	± 1 °C (according to DIN, EN, ISO and EPA methods)
Warm-up time:	from 20 °C to 160 °C within 10 minutes
Heating times:	4 pre-programmed heating times: 30 min, 60 min, 120 min, cont 8 free memory locations for individual time settings
Time range:	0h:01min – 9h:59min (1 min increments)
Safety:	Replaceable safety covers as contact protection Lockable protective lids Overheating protection
Interface:	Bidirectional serial RS 232, USB A (master) and USB B (slave) Allows connection of <b>NANOCOLOR®</b> T-Set (REF 919 917) (optional fully automatic calibration and preparation of test certificates for instrument control and monitoring)
Update:	Via Internet and USB stick
Power supply:	110 – 230 V ~, 50 / 60 Hz
Power consumption:	250 / 500 VA ( <b>NANOCOLOR® VARIO 4</b> ) or 125 / 250 VA ( <b>NANOCOLOR® VARIO C2</b> )
Dimensions:	290 x 287 x 146 mm ( <b>NANOCOLOR® VARIO 4</b> ) or 169 x 282 x 146 mm ( <b>NANOCOLOR® VARIO C2</b> )
Weight:	approx. 3.2 kg ( <b>NANOCOLOR® VARIO 4</b> ) or approx. 2.0 kg ( <b>NANOCOLOR® VARIO C2</b> )
Marking:	CE
Warranty:	2 years

### Declaration of conformity:

 These instruments conform to the following directives:  
2006/95/EC – Low-Voltage Directive  
2004/108/EC – EMV Directive

## 2. Operation

### 2.1. Initial operation

Remove the heating block and all accessories from the packing case. Place the heating block at room temperature on a dry and level surface, protected from direct sunlight. Connect the heating block to the mains supply and switch on the instrument with the main switch on the back.


#### Info

Serial number, actual software version and bootloader version are displayed. Press the info screen to return to the start screen.

### 2.2. Touchscreen-Display

The heating block is equipped with an illuminated coloured touchscreen display. The display is used by touching specific buttons. Familiarize yourself with the functions of the touchscreen by pressing individual buttons with your finger or a special pen. After touching a button, the relevant function is activated. Do not use sharp or pointed instruments to operate the touchscreen. Clean the touchscreen display with a micro-fibre cloth.

### 2.3. Program selection

Selection of a pre-programmed temperature. The temperature used last is always displayed first. Press the buttons repeatedly to scroll through the temperatures. Only programmed temperatures are displayed. The temperature values are always displayed in ascending order.

Selection of a pre-programmed heating time. The time used last is always displayed first. Press the buttons repeatedly to scroll through the heating times. Only programmed times are displayed. The time values are always displayed in ascending order.


Start of a selected program.

While warming up, the selected temperature is flashing in the display.

When the desired temperature is reached, a signal tone sounds and the selected time starts to run.

While warming up, the current temperature can be displayed by pressing the  button.

A running program can be stopped at any time. A program stop is triggered by pressing the  button.


*STOP* appears in the display.

Stop must be confirmed within 5 s by pressing the  button again, otherwise the program continues to run.

At the end of a program, a signal tone sounds and *END* appears in the display.

Press one of the following buttons (, , ) to call up the last program again.

## Graphic-mode


Press the  to zoom in and out.

Press the  button to go back to the start display.

If you use the **NANOCOLOR® VARIO 4** in the graphic mode (see menu, graphic-mode), press  to switch between the temperature curves of the two blocks.

### 3. Menu

General buttons:

**Esc**

Back to the main menu without confirming or saving changes.

**OK**

Changes are confirmed and/or saved.

**Language**


Select a language and confirm with .

**Graphic-mode**

A window is displayed to enable  or disable  the graphic mode. Select the required setting and confirm with .

## Temperature programming

The screenshot shows a temperature programming interface. At the top left, there is a 'P' icon and a temperature icon with an upward arrow. Below this is a grid of temperature slots. The first row contains 70°, 100°, 120°, and 148°, each with a small lock icon. The second row contains 160° and three empty slots with ellipses. The third row contains four empty slots with ellipses. At the bottom, there is a status bar with the date and time '01-05-2010 12:15:01' and a speaker icon.

Select an empty slot by pressing .

The screenshot shows a numeric keypad for temperature programming. At the top left, there is a 'P' icon and a temperature icon with an upward arrow. Below this is a grid of numeric keys. The first row contains 1, 2, and 3. The second row contains °, 4, 5, and 6. The third row contains 7, 8, and 9. The fourth row contains C, 0, and a checkmark icon. At the bottom, there is a status bar with the date and time '01-05-2010 12:15:01' and a speaker icon.


Enter the desired temperature in the range from 40 to 160 °C and confirm with .

User defined temperatures can be deleted or edited in the same way.

Pre-programmed temperatures can not be deleted or edited.

## Heating time programming

Select an empty slot by pressing .

Enter the desired heating time in the range from 0h:01min – 9h:59min and confirm with .

User defined heating times can be deleted or edited in the same way.

Pre-programmed heating times can not be deleted or edited.

## Date / time


Enter the time in the format hh:mm. Select the field for the date by touching. Enter the date in the format dd-mm-20yy and confirm with .

### Protective lid warning

A window is displayed to enable  or disable  the protective lid warning. Select the required setting and confirm with .

### Display

Press   to adjust brightness , contrast  and saturation  to the surrounding light conditions.


### Sound

Press   to adjust the sound level. The icon  in the display shows the adjusted sound level.

### T-Set

Press  to select the calibration mode. The calibration mode allows you to calibrate all programmed temperatures. After calibration, press  to confirm and save the new calibration data.

Press  to select the test mode. The test mode allows you to test all programmed temperatures, the calibration is not changed.

Press  to restore factory settings.

If you use the *NANOCOLOR® VARIO 4*, press   to select one block for calibration or testing.

Press  to start the program.

Deviations in temperature are displayed in different colors:

**Deviation**                      **Color**

< 1 °C                              green

1 < x < 2 °C                      yellow

> 2 °C                                red

Press  to save the calibration or test protocol on a USB stick.

**System settings** 

A window is displayed to restore the following factory settings.

User-defined temperatures

User-defined heating times

Display settings

Sound settings

Select the required setting and confirm with .

**Error list** 

A window with an error list is displayed.

Disposal in accordance with EU Directive 2002/96/CE

In compliance with local and national legal regulations (EU Directive 2002/96/CE), the MACHEREY-NAGEL company disposes old instruments free of charge.

Note: With effect from August 2005, disposal using public waste disposal facilities is no longer permitted. In the case of disposal, please contact your MACHEREY-NAGEL representative.

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# MACHERY-NAGEL

Your competent partner in analytical chemistry

## Tradition and Modernity – 100 years of experience

- Worldwide operating German company (founded 1911)
- Subsidiaries in France, Switzerland and USA
- Distributors in over 150 countries
- Longtime tradition in filter papers
- Development, production and sales of special products for water, environmental and food analysis, for biotechnology, chemical and pharmaceutical industry and medical diagnostics

## Quality and Diversity – 5 product ranges with over 25 000 products „Made in Germany“

### Filtration

- Filter Papers
- Filter Membranes
- Extraction Thimbles

### Rapid Tests

- Test Papers and Test Strips
- Urine Test Strips

### Water Analysis

- Colorimetric and Titrimetric Test Kits
- Photometric Water Analysis
- Microbiology

### Chromatography

- High Performance Liquid Chromatography (HPLC)
- Gas Chromatography (GC)
- Thin Layer Chromatography (TLC)
- Sample Preparation (SPE)

### Bioanalysis

- Kits for Purification of Nucleic Acids
- Kits for Purification of Proteins
- Transfer Membranes

## Service

- Technical support by our Customer Service Center
- Product-specific quality certificates and Material Safety Data Sheets
- Customised special products
- Easy ordering and short delivery times
- Large number of catalogues, flyer and product information – of course free of charge

## Contact

Technical Support and Customer Service Center for Water Analysis, Filtration, Rapid Tests, Medi-Test:

+49 (0)-2421 / 969 138

+49 (0)-2421 / 969 161

+49 (0)-2421 / 969 174

+49 (0)-2421 / 969 187

+49 (0)-2421 / 969 168

+49 (0)-2421 / 969 166

Product Management Water Analysis:

Product Management Filtration and Rapid Tests:

Please visit our website:

[www.mn-net.com](http://www.mn-net.com)

[www.mn-net.com](http://www.mn-net.com)

# MACHERY-NAGEL

MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren

Germany and international:

Tel.: +49 (0) 24 21 96 90

Fax: +49 (0) 24 21 96 91 99

e-mail: [sales-de@mn-net.com](mailto:sales-de@mn-net.com)

Switzerland:

Tel.: +41 (0) 62 388 55 00

Fax: +41 (0) 62 388 55 05

e-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

France:

Tel.: +33 (0) 3 88 68 22 68

Fax: +33 (0) 3 88 51 76 88

e-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)